


Assay values AC-Diff Controls

LOT	:	150
	:	2010, December 19th
REF	:	3892, 3893, 3894

Assay values and expected ranges for: Swelab AC900 series


	3892		3893		3894	
US UNITS	LOW		NORMAL		HIGH	
Parameters	Mean	Range ±	Mean	Range ±	Mean	Range ±
RBC 10 ¹² /l	2.28	0.12	3.61	0.15	5.29	0.18
Hb g/dl	5.8	0.3	11.3	0.4	17.9	0.5
Hct %	16.0	1.5	30.0	2.0	48.1	2.5
MCV fl	70	4	83	4	91	4
MCH pg	25.4	2.0	31.3	2.0	33.8	2.0
MCHC g/dl	36.3	2.0	37.7	2.0	37.2	2.5
RDW-CV %	19.2	2.0	15.9	2.0	16.2	2.0
Plts 10 ⁹ /l	65	15	170	20	349	40
WBC 10 ⁹ /l	2.8	0.4	8.7	0.6	23.9	1.8
NEUT%	69.8	5.0	50.4	4.0	21.2	5.0
LYMP %	17.0	6.0	37.4	3.0	69.6	3.0
MID %	13.2	3.0	12.2	3.0	9.2	3.0
NEUT#	2.0	0.2	4.4	0.4	5.1	1.2
LYMP #	0.5	0.3	3.3	0.3	16.6	0.7
MID #	0.4	0.1	1.1	0.2	2.2	0.7

Assay values and expected ranges for: Swelab AC900 series

	3892		3893		3894	
SI UNITS	LOW		NORMAL		HIGH	
Parameters	Mean	Range ±	Mean	Range ±	Mean	Range ±
RBC 10 ¹² /l	2.28	0.12	3.61	0.15	5.29	0.18
Hb mmol/l	3.6	0.2	7.0	0.3	11.1	0.3
Hct l/l	0.160	0.015	0.300	0.020	0.481	0.025
MCV fl	70	4	83	4	91	4
MCH fmol	1.58	0.12	1.94	0.12	2.10	0.12
MCHC mmol/l	22.6	1.2	23.4	1.4	23.1	1.6
RDW-CV %	19.2	2.0	15.9	2.0	16.2	2.0
Plts 10 ⁹ /l	65	15	170	20	349	40
WBC 10 ⁹ /l	2.8	0.4	8.7	0.6	23.9	1.8
NEUT%	69.8	5.0	50.4	4.0	21.2	5.0
LYMP %	17.0	6.0	37.4	3.0	69.6	3.0
MID %	13.2	3.0	12.2	3.0	9.2	3.0
NEUT#	2.0	0.2	4.4	0.4	5.1	1.2
LYMP #	0.5	0.3	3.3	0.3	16.6	0.7
MID #	0.4	0.1	1.1	0.2	2.2	0.7



Mallinckrodt Baker B.V.
P.O. Box 1 - 7400AA Deventer - Holland
Tel.: +31 (0)570 687500



The devices as mentioned in this insert comply with the In Vitro Diagnostic Medical Devices Directive 98/79/EC.

Mallinckrodt Baker B.V.
Teugseweg 20
P.O.Box 1
7400 AA Deventer
Holland
t: +31 570 687500
f: +31 570 687574
e-mail: jtbaker.nl@covidien.com

Mallinckrodt Baker, UK
107/112 Leadenhall Street
London
EC3A 4AH
United Kingdom
t: +44 1908 506000
f: +44 1908 503290
e-mail: jtbaker.uk@covidien.com

Mallinckrodt Baker, Italy
Via P.Portaluppi 11/2
20138 Milano
Italy
t: +39 02 580891
f: +39 02 58019264
e-mail: jtbaker.it@covidien.com

Mallinckrodt Baker, Germany
Im Leuschnerpark 4
D-64347 Griesheim
Germany
t: +49 6155 667939
f: +49 6155 667940
e-mail: jtbaker.de@covidien.com

Mallinckrodt Baker, France
33 Avenue Philippe Auguste
75011 Paris
France
t: +33 144 64 89 71/2
f: +33 144 64 89 73
e-mail: jtbaker.fr@covidien.com